

Report # _____

Company: NuMega
 Address: 11526 Sorrento Valley Road Ste. B2
San Diego, CA 92121
 Telephone: (858)793-6057 Fax: (858)796-2607
 Name: John Doe Date: Feb 2, 2017
 PO# _____ Email: NuMegalab@numegalabs.com

NuMega
 11526 Sorrento Valley Road Ste.B2
 San Diego, CA 92121
 Tel: (858)793-6057 Fax: (858)793-2607
<http://www.Numegalabs.com>
 Numegalab@numegalabs.com

PDF and FID

Sample ID		NMR [†]		Mass [‡]		Elemental Analysis Theoretical % please continue *			
		H1,P31,C13 etc.	Solvent	Range to scan	M.W.	C%	H%	N%	S%
A	A101	H1, P31	DMSO						
B	A102	H1, C13	CDCl3						
C									
D									
E									
F									
G									
H									

* **Elemental Analysis** (2mg per run) **Single:** **Duplicate:** **Or Duplicate if not in range:** _____ %

[†]**NMR (PDF):** -0.5 to 10 ppm, integration, peak picking for H-1, 0-200 ppm, peak picking for C-13. (**FID:** please request) Chemical shift (H-1). TMS (0), DMSO (2.50), CDCl₃(7.26), or D₂O (4.78).

[‡]**Mass:** 1mg in 1mL or 20-50 μM , positive and negative modes. Please specify MW range.

Your Instructions /Remarks: (Hygroscopic, Explosive, Volatile, Toxicity, Storage, EA theoretical %, etc.)

Confidentiality: NuMega maintains spectral information in strict confidence

Comment from Analyst: _____ **For NuMega Use Only**

Routine NMR	__ x \$45 =	
Non-Routine	__ x __ =	
Others	__ x __ =	

MASS	__ x \$35 =	

CHN (S/D)	__ x \$27/42 =	
S (S/D)	__ x \$27/42 =	

Analyst: _____ Date _____ Total Charge: \$ _____